

# Polymorph Characterization Using the Thermo Scientific Nicolet iS50 Raman Module

## Key Words

Crystallinity by Raman, Drug Discovery, Pharmaceuticals, Polymorphism, Raman Spectroscopy

Solid dose drugs comprise a considerable segment of the pharmaceutical market. The solid dose form can offer long shelf life, convenient drug delivery, and can be packaged into capsules, tablets, gel caps, and patches. *Polymorphism* is the ability of a molecule to exist in more than one stable or meta-stable crystalline state. Distinct crystalline forms commonly exhibit differences in pharmacologically relevant properties of a single active pharmaceutical ingredient (API) in its solid dose form. These properties include bioavailability, solubility, rate of absorption, duration, stability, manufacturability, and processability. Manufacturing factors also can influence polymorphism and can include solvents used in product formulation, the rate of product recrystallization, salt selection, and temperature. Since polymorphism can have profound effects on an API's pharmacodynamics and manufacturing process, it becomes critical to fully characterize and understand polymorphic behavior.

The information-rich content of Raman spectra and the availability of reliable commercial Raman instruments have given rise to the use of Raman spectroscopy in many aspects of drug development<sup>1-3</sup>. Raman spectroscopy has proven to be a valuable analytical method for differentiating between different polymorphs, salts, co-crystals, solvates, and hydrates of potential drug candidates. In many cases, the Raman spectra can distinguish between many of these existing forms as well as potentially identifying new, unknown polymorphs. Raman analysis offers several benefits including the lack of special sample preparation for solid dosage forms and the ability to sample through packaging, like blister packs, without interference from the clear packaging material.



In this application note we describe using the Thermo Scientific<sup>™</sup> Nicolet<sup>™</sup> iS<sup>™</sup>50 Raman module in the Thermo Scientific Nicolet iS50 FT-IR spectrometer sampling compartment to characterize different polymorphs of a model compound, 5-Methyl-2-[(2-nitrophenyl)amino]-3-thiophenecarbonitrile (ROY)<sup>4-6</sup>. This compound was kindly provided by Professor Lian Yu, of the University of Wisconsin, School of Pharmacy, Division of Pharmaceutical Sciences. This unique compound exhibits different colors for different polymorphs, and we describe here the analysis of three ROY polymorphs (Red, Orange, and Yellow) seen in Figure 1.




Figure 1: Pictures of the various polymorphs found for 5-Methyl-2-[(2-nitrophenyl)amino]-3-thiophenecarbonitrile (ROY)

### Experimental Methods

The spectra described in this paper were acquired with the Nicolet iS50 Raman Module (Thermo Fisher Scientific, Madison WI). The Raman module mounts in the sample compartment of the Nicolet iS50 spectrometer and is equipped with a 1064 nm near-infrared laser which significantly reduces fluorescence interference relative to visible laser excitation. The samples were mounted in a mini-well plate that was placed on the fully automated X, Y, Z stage of the Raman module. Figure 2 shows a combined image of the wells containing the three samples that was acquired with the digital camera embedded in the module. This provided a video record of the precise position on the sample where the spectrum was acquired. The total data collection and analysis time for samples such as this depends upon the method employed. Using a single “aim-focus-collect” process (convenient for 3 samples) required about 3 minutes. The Array Automation software used here automates the analysis for a grid of samples, permitting unattended operation. Array Automation also links to spectral analysis tools, allowing the experiment to move from data collection to final report automatically.




Figure 2: This composite image from the iS50 Raman module shows the three ROY polymorphs in three 1 mm diameter wells in a mini-plate. Starting from the left; Red, Yellow and Orange polymorphs




Figure 3: Full Raman spectra from three polymorphs of ROY




Figure 4: Expansion of the nitrile peak showing shift for three different ROY polymorphs

The spectra from the three samples contain a number of significant spectral differences that are shown in Figure 3. One clear difference is the position of the nitrile peak near 2200 cm<sup>-1</sup> that has been expanded in Figure 4. This peak shift appears to correlate with the different angle observed for the nitrile functional group in the different crystal forms.

Prolonged exposure of the Orange polymorph to the excitation laser produced an interesting effect. Spectra collected over time at high laser power showed an increased presence of a second nitrile peak that closely matched the peak from the Yellow polymorph. Figure 5 shows the resulting mixture spectra on the bottom with the three pure polymorphs displayed on top. We applied the Thermo Scientific OMNIC™ curve analysis software to fit the mixture spectra. The results of one curve analysis are shown in Figure 6, demonstrating that Raman spectra can be used to determine the amount of different polymorphs present in a sample. The Thermo Scientific OMNIC TQ Analyst™ multivariate statistical analysis software would enable further analysis of this sample, including quantitative or qualitative analyses.

Raman spectroscopy provides pharmaceutical companies with a versatile tool for identifying the proper polymorph (or polymorph blend) during drug discovery. This allows optimization of the production process to yield the desired product and speeds scale-up work into manufacturing. Post-launch, Raman assists in identifying counterfeit drugs or patent infringement issues to protect the intellectual property involved. Beyond pharmaceuticals, many laboratories – from research labs to forensic labs – can benefit from the additional information Raman spectroscopy offers in understanding the properties of solid dose drugs.




Figure 5: A comparison of the spectra from the pure materials (top) with spectra from samples containing a mixture of two polymorphs




Figure 6: Results of curve fitting that provide a good estimate of relative concentrations from peak heights or areas for each polymorph

## Conclusion

Optimizing the choice of polymorphs of an API has become a significant step in the drug development process. The polymorphism and salt form of an API greatly affect a drug's performance and stability, drastically changing absorption and other characteristics *in vivo*. During the discovery process, identifying and comparing the properties of different polymorphs, salts, solvates, complexes, and co-crystals becomes vital for building intellectual property and creating a competitive advantage for pharmaceutical companies.

The Nicolet iS50 FT-IR spectrometer equipped with the iS50 Raman module (Figure 7) can be extremely valuable for characterizing polymorphs and mixtures during drug development as well as providing a tool for identifying counterfeit drugs or patent infringement. The fully integrated system – with Raman and FT-IR, plus other optional tools – offers an ideal workstation for drug development, including screening polymorphic forms. Raman analysis is fast, requires little to no sample preparation and can be used to characterize bulk or single crystal samples. Raman permits measurements through many packaging materials such as blister packs. The iS50 Raman module enables experiment automation, and does not require consumables or extra safety precautions. Raman spectroscopy using the iS50 Raman module represents an economical, reliable and complementary technique to DSC (differential scanning calorimetry) or X-ray diffraction for characterizing APIs.



Figure 7: The Thermo Scientific Nicolet iS50 spectrometer with the iS50 Raman module in the sample compartment with a well plate designed for analyzing solid dose tablets

## References

- Lowry, S.; Dalrymple, D.; Song, A.; Rosso, V.; Pommier, C.; Venit, J. Integrating a Raman Microscope into the Workflow of a High-Throughput Crystallization Laboratory. *J. Assoc. Lab. Autom.* 2006, 11(2), 75–84.
- Igo, D.; Chen, P. Vibrational Spectroscopy of Solid-State Forms – Application and Examples. In *Applications of Vibrational Spectroscopy in Pharmaceutical Research and Development*; Pivonka, D., Chalmers, J., Griffiths, P. Eds; John Wiley & Sons: New York, 2007; pp 293–308.
- Bugay, D.; Henck, J.; Longmire, M.; Thorley, F. Raman Analysis of Pharmaceuticals. In *Applications of Vibrational Spectroscopy in Pharmaceutical Research and Development*; Pivonka, D., Chalmers, J., Griffiths, P. Eds; John Wiley & Sons: New York, 2007, pp 239–262.
- Yu, L. Red, Orange and Yellow Polymorphs of a Conformationally Flexible Molecule: Single-Crystal Spectroscopy, Optical Crystallography and Computational Chemistry. *J. Phys. Chem. A* 2002, 106, 544–550.
- Chen, S.; Guzei, I. A.; Yu, L. New Polymorphs of ROY and New Record for Coexisting Polymorphs of Solved Structures. *J. Am. Chem. Soc.* 2005, 127, 9881–9885.
- Yu, L. Polymorphism in Molecular Solids: An Extraordinary System of Red, Orange, and Yellow Crystals. *Acc. Chem. Res.* 2010, 43, 1257–1266. Cover article.

## www.thermoscientific.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

<b>Africa</b> +27 11 822 4120	<b>Denmark</b> +45 70 23 62 60	<b>India</b> +91 22 6742 9434	<b>New Zealand</b> +64 9 980 6700
<b>Australia</b> +61 3 9757 4300	<b>Europe-Other</b> +43 1 333 50 34 0	<b>Italy</b> +39 02 950 591	<b>Russia/CIS</b> +43 1 333 50 34 0
<b>Austria</b> +43 1 333 50 34 0	<b>Finland/Norway/Sweden</b> +46 8 556 468 00	<b>Japan</b> +81 45 453 9100	<b>Spain</b> +34 914 845 965
<b>Belgium</b> +32 53 73 42 41	<b>France</b> +33 1 60 92 48 00	<b>Latin America</b> +1 561 688 8700	<b>Switzerland</b> +41 61 716 77 00
<b>Canada</b> +1 800 530 8447	<b>Germany</b> +49 6103 408 1014	<b>Middle East</b> +43 1 333 50 34 0	<b>UK</b> +44 1442 233555
<b>China</b> +86 10 8419 3588		<b>Netherlands</b> +31 76 579 55 55	<b>USA</b> +1 800 532 4752

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific