

Elemental analysis: CHNS characterization of rocks by flash combustion

Author: Liliana Krotz, Thermo Fisher Scientific, Milan, Italy

Keywords: CHNS, combustion, NCS, rocks, TOC

Goal

To assess the performance of the elemental analyzer for CHNS determination of rocks in terms of the accuracy, precision, and repeatability.

Introduction

Elemental analysis is used for the characterization of rocks. The determination of Nitrogen and Carbon in rocks is important for the evaluation of organic matter. The differentiation of Total Carbon and Total Organic Carbon as the determination of sulfur are also useful parameters to characterize rocks.

Traditional methods are no longer suitable for routine analysis of rocks, due to their time-consuming preparation and the use of environmentally hazardous reagents. For these reasons, the use of an accurate instrumental analytical techniques is required. As the demand for improved sample throughput, reduction of operational costs and minimization of human errors is increasing notably, a simple and automated technique, which allows fast analysis with an excellent reproducibility is the key for the elemental determination of rocks.



The Thermo Scientific™ FlashSmart™ Elemental Analyzer (Figure 1), based on the dynamic flash combustion of the sample, copes effortlessly with the wide array of laboratory requirements such as accuracy, day to day reproducibility and high sample throughput. The FlashSmart EA allows the automated elemental determination of rock samples and the same analytical conditions can be used for the differentiation between the Total Carbon and Total Organic Carbon determination after an acid pre-treatment of the sample. Through its flexibility, the FlashSmart EA allows also the simultaneous NCS analysis while for trace Sulfur determination, the analyzer has been coupled with the flame photometric detector (FPD). This method combines the advantages of the elemental analyzer with the sensitivity, selectivity, and robustness of the FPD Detector.



Figure 1. Thermo Scientific™ FlashSmart™ Elemental Analyzer.

Methods

For CHNS abundance determination, the FlashSmart EA operates with the dynamic flash combustion of the sample. Samples are weighed in tin containers and introduced into the combustion reactor via the Thermo Scientific™ MAS Plus Autosampler alongside a pulse of oxygen. After combustion, the produced gases are carried in a helium carrier gas to a layer filled with copper. The analyte then enters the GC column, which separates the produced gases before detection by a Thermal Conductivity Detector (TCD) (Figure 2). For weight percent determination a complete report is automatically generated by the Thermo Scientific™ EagerSmart™ Data Handling Software and displayed at the end of the analysis. For S (single determination) or simultaneous NCS configuration, after combustion of the sample the resultant gases are carried by a helium flow to a layer filled with copper, then through a water trap, a GC column and finally, detected by the thermal conductivity detector (TCD) (Figure 3), while for trace sulfur analysis, after the water trap, the gases are carried by a helium flow through a short GC column and finally, detected by the Flame Photometric Detector (FPD), (Figure 4).




Figure 2. CHNS configuration.




Figure 3. Single Sulfur or NCS configuration.




Figure 4. Sulfur configuration by FPD Detector.

The differentiation of Total Carbon (TC) and Total Organic Carbon (TOC) was performed by sample manipulation prior to analysis following the Official Italian Method on Soils Analytical Chemistry, Method 248 (Gazzetta Ufficiale).

TOC was determined after removing carbonates by acidification of the sample with HCl 1:1 (Figure 5) using the kit showed in Figure 6. The two analyses TC and TOC were performed consecutively using the same analytical conditions of the instrument.




Figure 5. Method for TOC determination.



Figure 6. Kit for Total Organic Carbon (TOC) for solid samples.

For CHNS abundance determination, the calibration curve was produced by analyzing 2–3 mg BBOT and using the K factor as the calibration method. The rock samples were analyzed 10 times to evaluate the repeatability. The trace sulfur content of sample code 11 and 14 was determined using the FlashSmart EA coupled with the Flame Photometric Detector. Table 1 shows the sample information (the thermal maturity increases approximately with age in these samples), and the sample weight used for CHNS and TOC determination. Table 2 shows the CHNS and TOC data with the relative RSD% obtained for each sample.

Results

The analysis of 15 rock samples with different geological ages and thermal maturity, including 8 USGS* rock standards, were performed to demonstrate the performance of the Analyzer. Samples were homogenized by a ball mill.

*USGS: the United States Geological Survey is a government organization that studies the geological history of the United States and provides analytical reference materials.

Table 1. Rock sample information.

| Code | Rock sample name | Geological unit origin | Age (billion of years) | Weight (mg) for CHNS | Weight (mg) for TOC |
|------|------------------|-----------------------------------|------------------------|----------------------|---------------------|
| 1 | USGS SGR-1 | Green River Shale, USA | 0.05 | 3–4.5 | 3–4.5 |
| 2 | USGS SDO-1 | Devonian Ohio Shale, USA | 0.37 | 7–10 | 8–10 |
| 3 | USGS SHWFD-1 | Woodford Shale, USA | 0.36 | 7–10 | 8–10 |
| 4 | USGS SHBOQ-1 | Boquillas Shale, USA | 0.07 | 7–8 | 8–10 |
| 5 | USGS SCO-1 | Cody Shale, USA | 0.07 | 15–20 | 8–10 |
| 6 | USGS SBC-1 | Brush Creek Shale, USA | 0.31 | 15–25 | 8–10 |
| 7 | USGS BHVO-2 | Hawaiian Basalt, USA | Modern | 15–25 | 8–10 |
| 8 | USGS SDC-1 | Mica Schist, USA | Unknown | 95–105 | 8–10 |
| 9 | MR21011 | Mt McRae Shale, Australia | 2.50 | 10–15 | 8–10 |
| 10 | J18 | Jeerinah Formation, Australia | 2.66 | 10–15 | 8–10 |
| 11 | NS1282-1 | Nonesuch Shale, USA | 1.10 | 15–25 | 8–10 |
| 12 | SC20-1 | Sheep Creek, Belt Supergroup, USA | 1.45 | 15–20 | 8–10 |
| 13 | SC20-51 | Sheep Creek, Belt Supergroup, USA | 1.45 | 15–25 | 8–10 |
| 14 | MC-4 | Mosquito Creek Group, Australia | 2.85 | 15–25 | 8–10 |
| 15 | T37 | Tumbiana Formation, USA | 2.72 | 15–25 | 8–10 |

Table 2. CHNS and TOC data of rocks.

| Code | CHNS abundance and RSD% | | | | | | | | TOC abundance and RSD% | |
|------|-------------------------|------|--------|------|--------|------|--------|------|------------------------|------|
| | N% | RSD% | C% | RSD% | H% | RSD% | S% | RSD% | TOC% | RSD% |
| 1 | 0.883 | 0.77 | 27.84 | 0.37 | 3.17 | 0.52 | 1.47 | 1.22 | 24.56 | 0.25 |
| 2 | 0.357 | 0.41 | 9.62 | 0.33 | 1.46 | 1.25 | 5.29 | 0.27 | 9.33 | 0.09 |
| 3 | 0.275 | 0.81 | 7.93 | 0.49 | 0.944 | 0.95 | 1.03 | 0.82 | 7.70 | 0.49 |
| 4 | 0.102 | 1.99 | 11.55 | 0.90 | 0.747 | 1.28 | 1.57 | 1.11 | 5.27 | 0.55 |
| 5 | 0.0492 | 1.28 | 1.06 | 1.04 | 0.608 | 1.04 | 0.0219 | 1.02 | 0.2831 | 1.15 |
| 6 | 0.0578 | 0.75 | 2.08 | 0.85 | 0.785 | 0.78 | 0.259 | 0.83 | 1.11 | 0.40 |
| 7 | 0.0008 | 6.04 | 0.0221 | 2.96 | 0.0246 | 2.75 | 0.0092 | 2.24 | 0.0198 | 1.98 |
| 8 | 0.0030 | 3.08 | 0.0547 | 1.78 | 0.219 | 1.56 | 0.0551 | 1.24 | 0.0331 | 1.85 |
| 9 | 0.0966 | 0.97 | 5.95 | 0.87 | 0.878 | 0.86 | 9.76 | 0.63 | 5.62 | 0.33 |
| 10 | 0.0114 | 1.87 | 3.06 | 0.72 | 0.869 | 0.78 | 2.65 | 0.64 | 2.93 | 0.66 |
| 11 | 0.0418 | 1.71 | 0.0587 | 1.47 | 0.458 | 1.26 | 0.0033 | 4.25 | 0.0562 | 0.98 |
| 12 | 0.0251 | 1.98 | 0.584 | 0.51 | 0.272 | 1.44 | 0.342 | 0.96 | 0.5485 | 0.35 |
| 13 | 0.0069 | 2.74 | 0.430 | 0.90 | 0.756 | 1.34 | 0.114 | 1.17 | 0.4000 | 1.02 |
| 14 | 0.0030 | 3.54 | 0.179 | 0.71 | 0.654 | 0.87 | 0.0039 | 3.35 | 0.1482 | 1.05 |
| 15 | 0 | 0 | 0.470 | 0.30 | 0.336 | 0.41 | 0.379 | 0.33 | 0.3768 | 0.88 |

Other four rock samples were analyzed for NCS and TOC determination. The calibration curve was performed analyzing 2–3 mg BBOT and using the K factor as the calibration method. The rock samples were homogenized by a ball mill and analyzed in triplicate to evaluate the repeatability. For NCS determination, the sample weight was 5–10 mg for samples A and B, and 15–20 mg for samples C and D, while for TOC determination the sample weight was 5–6 mg. Table 3 shows the data obtained.

Table 3. NCS and TOC determination of rocks.

| Sample | N% | RSD% | C% | RSD% | S% | RSD% | TOC% | RSD% |
|--------|--------|------|---------|-------|--------|------|--------|------|
| Rock A | 0.0546 | 0.28 | 12.2622 | 0.001 | 1.0097 | 0.64 | 1.6601 | 0.48 |
| | 0.0544 | | 12.2625 | | 1.0095 | | 1.6543 | |
| | 0.0543 | | 12.2624 | | 1.0209 | | 1.6702 | |
| Rock B | 0.0540 | 1.17 | 5.2900 | 0.44 | 4.1788 | 0.57 | 4.6462 | 0.53 |
| | 0.0552 | | 5.3362 | | 4.1469 | | 4.6003 | |
| | 0.0550 | | 5.3191 | | 4.1327 | | 4.6381 | |
| Rock C | 0.0062 | 2.40 | 0.1791 | 2.05 | 0.1643 | 3.89 | 0.1511 | 2.20 |
| | 0.0065 | | 0.1733 | | 0.1774 | | 0.1559 | |
| | 0.0064 | | 0.1800 | | 0.1690 | | 0.1577 | |
| Rock D | 0.0033 | 6.06 | 0.0997 | 1.14 | 0.3677 | 0.10 | 0.0269 | 2.68 |
| | 0.0035 | | 0.0982 | | 0.3684 | | 0.0255 | |
| | 0.0031 | | 0.0975 | | 0.3682 | | 0.0260 | |

Three silicate rock samples were analyzed using CHNS configuration. The calibration curve was performed analyzing 2–3 mg BBOT and using the K factor as the calibration method. The samples were homogenized by a ball mill and analyzed in triplicate to evaluate the repeatability, sample weight 15–20 mg. Table 4 shows the data obtained.

Table 4. CHNS data of silica rock samples.

| Sample | N% | RSD% | C% | RSD% | H% | RSD% | S% | RSD% |
|--------|--------|------|--------|------|-------|------|--------|------|
| 1 | - | - | 0.0491 | 1.40 | 0.272 | 0.63 | 0.0082 | 1.90 |
| | | | 0.0498 | | 0.275 | | 0.0079 | |
| | | | 0.0478 | | 0.272 | | 0.0080 | |
| 2 | 0.0024 | 7.91 | 0.0136 | 3.03 | 0.128 | 1.66 | - | - |
| | 0.0027 | | 0.0128 | | 0.124 | | | |
| | 0.0028 | | 0.0132 | | 0.125 | | | |
| 3 | - | - | 0.0643 | 0.76 | 0.184 | 0.84 | - | - |
| | | | 0.0644 | | 0.181 | | | |
| | | | 0.0652 | | 0.183 | | | |

Finally, a rock sample was analyzed at about 100 mg in CHNS configuration using the large tin container to weight the sample to demonstrate the repeatability of ten replicates. The calibration curve was performed analyzing 2–3 mg BBOT and using the K factor as the calibration method. Table 4 shows the data obtained.

Table 5. CHNS data a rock sample at about 100 mg sample weight.

| Weight (mg) | N% | C% | H% | S% |
|-----------------|--------|--------|--------|--------|
| 101.5 | 0.0031 | 0.0559 | 0.223 | 0.0551 |
| 102.5 | 0.0030 | 0.0551 | 0.220 | 0.0543 |
| 101.6 | 0.0031 | 0.0557 | 0.222 | 0.0549 |
| 101.2 | 0.0031 | 0.0558 | 0.223 | 0.0550 |
| 102.9 | 0.0030 | 0.0549 | 0.219 | 0.0541 |
| 106.0 | 0.0029 | 0.0545 | 0.218 | 0.0562 |
| 100.6 | 0.0029 | 0.0536 | 0.215 | 0.0553 |
| 100.2 | 0.0029 | 0.0538 | 0.216 | 0.0555 |
| 102.2 | 0.0029 | 0.0544 | 0.218 | 0.0561 |
| 103.5 | 0.0029 | 0.0531 | 0.213 | 0.0548 |
| Average | 0.0030 | 0.0547 | 0.219 | 0.0551 |
| Std.Dev. | 0.0001 | 0.0010 | 0.0034 | 0.0007 |
| RSD% | 2.08 | 1.78 | 1.56 | 1.24 |

Conclusions

For the quantitative determination of nitrogen, carbon, hydrogen and sulfur and TOC, the all-in-one FlashSmart EA is the optimal solution in geology in sample matrices with a wide range of concentrations spanning low to high amounts. Specifically, the FlashSmart EA demonstrates excellent repeatability, reproducibility, accuracy, and precision, as automation, speed of analysis and cost per analysis.

No memory effect was observed when changing the sample type, indicating the complete conversion and detection of all elements.

Thanks to the modularity of the FlashSmart EA, the hardware, autosamplers and software can be readily used for other configurations such as CHN/O, CHN/S, CHNS/O, CHNS/CHNS, CHN/CHN, NC (single reactor)/S, N-Protein (single reactor)/S and more. This can be achieved only by changing consumables.

The USGS rocks standards cover a wide range in CHNS abundances, making them an ideal suite of reference materials for geochemical studies.

Acknowledgments

We would like to thank Eva Stüeken, School of Earth and Environmental Sciences, St. Andrews University, Scotland, UK for providing the rock samples and collaboration in this project; Mattia Iammarino, IIS Curie-Sraffa student (Milano, Italy) for his collaboration in the analysis of the samples during his internship in Thermo Scientific OEA Demo Lab (Milan, Italy).

Find out more at thermofisher.com/OEA