

MOBILE RACKS

CAGES, RACKS FOR CONVENTIONAL HUSBANDRY

Take advantage of ZOONLAB's 50 years of experience in the manufacture of mobile racks for the storage of cages. Our selection and processing of section tubes, particularly with hairpin bracket design, all dirt traps and surfaces prone to dust accumulation will be avoided. Fitted with autoclavable castors with a diameter of 100 mm,

2 of them are equipped with a locking device, the racks can be positioned in any desired position. The racks are manufactured using premium stainless steel (AISI 304). The castor housing is also manufactured using stainless steel. We also offer autoclavable cover hoods as an optional extra for the racks.





RACKS WITH SWIVEL CASTORS



Rack for 30 type 2L cages, single-sided



Rack for 15 type 4 cages, single-sided


Rack for	Number of levels	Single-sided/double-sided	Dimensions (W x D x H) in mm	Order no.
Type 1 cage, long version – hairpin consoles				
42 type 1L cages	6 levels, each containing 7 cages	Single-sided	1,450 x 510 x 1,800	5033801
49 type 1L cages	7 levels, each containing 7 cages	Single-sided	1,450 x 510 x 1,800	5033800
84 type 1L cages	6 levels, each containing 14 cages	Double-sided	1,450 x 800 x 1,800	5033803
98 type 1 L cages	7 levels, each containing 14 cages	Double-sided	1,450 x 800 x 1,800	5033802
Type 2 cage – hairpin consoles				
30 type 2 cages	6 levels, each containing 5 cages	Single-sided	1,285 x 510 x 1,800	5033806
35 type 2 cages	7 levels, each containing 5 cages	Single-sided	1,285 x 510 x 1,800	5033807
60 type 2 cages	6 levels, each containing 10 cages	Double-sided	1,285 x 800 x 1,800	5033804
70 type 2 cages	7 levels, each containing 10 cages	Double-sided	1,285 x 800 x 1,800	5033805
Type 2 cage, long version – hairpin consoles				
30 type 2 L cages	6 levels, each containing 5 cages	Single-sided	1,285 x 510 x 1,800	5033810
35 type 2 L cages	7 levels, each containing 5 cages	Single-sided	1,285 x 510 x 1,800	5033811
60 type 2 L cages	6 levels, each containing 10 cages	Double-sided	1,285 x 800 x 1,800	5033809
70 type 2 L cages	7 levels, each containing 10 cages	Double-sided	1,285 x 800 x 1,800	5033808
Type 3 cage – hairpin consoles				
24 type 3 cages	6 levels, each containing 4 cages	Single-sided	1,285 x 510 x 1,800	5033814
28 type 3 cages	7 levels, each containing 4 cages	Single-sided	1,285 x 510 x 1,800	5033815
48 type 3 cages	6 levels, each containing 8 cages	Double-sided	1,285 x 800 x 1,800	5033813
56 type 3 cages	7 levels, each containing 8 cages	Double-sided	1,285 x 800 x 1,800	5033812
Type 3 H cage – hairpin consoles				
24 type 3 H cages	6 levels, each containing 4 cages	Single-sided	1,285 x 510 x 1,800	5033817
48 type 3 H cages	6 levels, each containing 8 cages	Double-sided	1,285 x 800 x 1,800	5033816
Type 4 cage – angle steel structure				
15 type 4 cages*	5 levels, each containing 3 cages	Single-sided	1,340 x 650 x 1,800	5033818
18 type 4 cages	6 levels, each containing 3 cages	Single-sided	1,340 x 650 x 1,800	5033819
20 type 4 cages*	5 levels, each containing 4 cages	Single-sided	1,760 x 650 x 1,800	5033820
24 type 4 cages	6 levels, each containing 4 cages	Single-sided	1,760 x 650 x 1,800	5033821

* also suitable for elevated lids
Subject to technical modifications.

TYPE 4 RACK WITH GRIP



Type 4 rack with 12 cages



Type	Description	Dimensions (W x D x H) in mm	Order no.
Type 4 rack	for 12 type 4 cages on 4 levels, made of stainless steel (AISI 304), manufactured in a welded design, stained and passivated in a full bath, equipped with 4 castors with a diameter of 100 mm complete with bolt hole, 2 of them are equipped with locking device, autoclavable. Clearance between the lower edge of the cage and the floor: 305 mm	1,224 x 654 x 1,797	5033585
Option			
Type 4 lid	complete with opening in order to insert a feed container, complete with a water bottle holder, made of stainless steel (AISI 304), electropolished		5023305
Feed container	for approx. 1.4kg of pellets		5083737



LECTERN

MADE OF STAINLESS STEEL, AUTOCLAVABLE

Type	Dimensions (W x D x H) in mm	Order no.
Base frame with document compartment	650 x 470 x 1,100	5115026
Options in combination with the base frame		
Printer compartment	591 x 447 x 300	5115027
Drawer	531 x 295 x 150	5115028
Extendable keyboard rack	531 x 300 x 20	5115029
Inclined rest for 12° inclination, applied loosely to the lectern	606 x 449 x 139	5115024

Subject to technical modifications.



VARIABLE RACK GV

In addition to the traditional rack types, Z00NLAB also supplies a variable rack which allows the accommodation of any desired cage type.

It is characterised by 6 extremely sturdy levels produced from grid elements (the grid bars run in the sliding direction of the cage) which are also infinitely heightadjustable as required.



The entire rack is made of stainless steel (AISI 304).

4 swivel castors with a diameter of 100 mm, 2 of them equipped with locking device. The castor housing is made of stainless steel and the castors are also autoclavable.

OPTION: Holder for feeding bowls
Further options available on request



Cage type	Capacity per level	Dimensions in mm (L x D x H)	Order no.
1 L	6	1,230 x 550 x 1,780	5033520
2	5		
2	10*		
HRC 500	6		
2 L	5		
3 / 3H	4		
4	3		
Option			
GV-holder for feeding bowls			5063398

* when placed on both sides

Subject to technical modifications.