

RABBIT CAGES

INDIVIDUAL CAGES FOR RABBITS

The individual ZONLAB cages for rabbits provide ideal conditions for animal husbandry in terms of size and shape and are also easy to use.

Depending on the material selected, the cage trays are heat-resistant up to max. 90°C or max. 120°C.





RABBIT CAGES

INDIVIDUAL CAGES FOR RABBITS

The ZOONLAB rabbit cages have been designed in accordance with the 2010/63/EU Directive of the European Parliament 22.09.2010 (formerly ETS 123).

They combine:

- A high level of flexibility through the interlockable racks that enable the floor space to be increased to twice the size or even more when the partition walls are removed (without the need for tools).
- The cage trays with a height of 10 cm as well as the animal dropping trays are available in the LURAN (heat-resistant up to max. 90°C) as well as in the autoclavable NORYL (heat-resistant up to max. 120°C) designs.
- The raised resting area (available in LURAN and NORYL designs) can be folded up which enables the cage tray to be removed easily without having to first remove the resting area.
- All metal parts are made of stainless steel (AISI 304).

- The rack castors have a diameter of 100 mm and two of them on each rack are equipped with locking device.
- The racks can optionally be equipped with an automatic nipple drinking device.



Racks and accessories for rabbits weighing up to 3 kg/habitat height of 50 cm or up to 5 kg/habitat height of 60 cm, floor space 4,200 cm².

Type	Material	Dimensions (W x D x H) in mm	Order no.
Rack GK42-100/50-3, 3 levels each containing 1 cage, 50 cm high incl. 3 grid doors	Stainless steel (AISI 304)	800 x 780 x 1,890	5053164
Rack GK42-100/60-2, 2 levels each containing 1 cage, 60 cm high incl. 2 grid doors	Stainless steel (AISI 304)	800 x 780 x 1,790	5053166
Rabbit cage K42/100-L	LURAN	730 x 730 x 100	3045338
Feed container, approx. 1,25 litres, FUKA	Stainless steel (AISI 304)	W = 125 mm	5083721
Bottle holder for 1 bottle 750 ml, FH 75/100-1	Stainless steel (AISI 304)	72 x 72 x 182	5121005
Water bottle, 750ml, with groove, F 75	Polycarbonate	72 x 72 x 180	3125004
Silicone sealing ring, SDR 1	Silicone	-	3125101
Bottle cap incl. 40 mm tube with ball closure, TK-40	Stainless steel (AISI 304)	-	3125042
Animal dropping tray, KK 4200-L	LURAN	750 x 730 x 40	3094229
Partition wall, non-transparent for 3-level rack, 42-50-T	TRESPA	755 x 430	4043762
Partition wall, non-transparent for 2-level rack, 42-60-T	TRESPA	755 x 530	4043758
Options for rabbit cages			
Raised resting area	NORYL	1,650 cm ²	5044930S
Raised resting area	LURAN	1,650 cm ²	5044931S
Rabbit cage K 42/100-N	NORYL	730 x 730 x 100	3045339
Water bottle, 750 ml, with groove, F 75 PSU	Polysulfone	72 x 72 x 180	3125016
Animal dropping tray, KK 4200-N	NORYL	750 x 730 x 40	3094230
Partition wall, transparent for 3-level rack, 42-50-P	Polycarbonate	755 x 430	on request
Partition wall, transparent for 2-level rack, 42-60-P	Polycarbonate	755 x 530	4043759
Automatic nipple drinking device for 3-level rack, AT 42/3	Stainless steel (AISI 304)	-	5054060
Automatic nipple drinking device for 2-level rack, AT 42/2	Stainless steel (AISI 304)	-	5054061

Subject to technical modifications.

RABBIT CAGES

INDIVIDUAL CAGES FOR RABBITS



GK54/2-60 – two interlocked cage racks

Cage racks and accessories for rabbits over 5 kg/habitat height of 60 cm, floor space 5,400 cm²

Type	Material	Dimensions (W x D x H) in mm	Order no.
Rack GK 54/2-60, 2 levels each containing 1 cage, 60 cm high incl. 2 grid doors	Stainless steel (AISI 304)	940 x 860 x 1,775	5053144
Rack GK 54/3-60, 3 levels each containing 1 cage, 60 cm high incl. 3 grid doors, top level can be removed	Stainless steel (AISI 304)	940 x 860 x 2,300	5053149
Rabbit cage K5400-N	NORYL	855 x 775 x 100	3045434
Feed container, approx. 1.25 litres, FUKA	Stainless steel (AISI 304)	W = 125 mm	5083721
Bottle holder for 1 bottle 750 ml, FH 75/100-1	Stainless steel (AISI 304)	–	5121005
Water bottle, 750 ml, with groove, F 75	Polycarbonate	72 x 72 x 180	3125004
Silicone sealing ring, SDR 1	Silicone	–	3125101
Bottle cap incl. 40 mm tube with ball closure, TK-40	Stainless steel (AISI 304)	–	3125042
Animal dropping tray, KK 5400-N	NORYL	855 x 775 x 40	3094228
Partition wall, non-transparent, 54-60-T	TRESPA	820 x 540	4043752
Options for rabbit cages			
Raised resting area	NORYL	1,945 cm ²	5044932S
Water bottle, 750 ml, with groove, F 75 PSU	Polysulfone	72 x 72 x 180	3125016
Partition wall, transparent, 54-60-P	Polycarbonate	820 x 540	4043753
Hay rack	Stainless steel (AISI 304)	–	5083789
Automatic nipple drinking device, AT 54/2	Stainless steel (AISI 304)	–	5054056

Subject to technical modifications.



RABBIT CAGES

INDIVIDUAL CAGES FOR RABBITS

As an expansion to our rabbit cage range, the rack offers you ideal handling when it comes to storage / cleaning or the transport of resting areas. The lightness and manoeuvrability of the rack meets all of the ergonomic requirements. The rack is made of stainless steel (AISI 304), is glass-bead blasted with autoclavable castors, two of them are equipped with locking device.



Type	Material	Dimensions (W x D x H) in mm	Order no.
Rack for the cleaning / storage and transportation of raised resting areas.	Stainless steel (AISI 304)	742 x 512 x 1,780	5033561

Subject to technical modifications.