

# DR1900 PORTABLE SPECTROPHOTOMETER

For Reliable Measurements Wherever You Are

Whatever you need to measure, the DR1900 Spectrophotometer from Hach provides accurate results. Wherever you need to measure, this robust, portable device is always to hand.



## Key benefits

- 200+ pre-programmed methods
- 50 freely programmable user applications
- Intuitive and user friendly
- Very robust in the field and in the laboratory
- Protection class IP67
- Compact dimensions and lightweight
- Available with optional outdoor backpack

Detailed information on technical data and accessories can be found on our website.

## Portable and reliable

The DR1900 combines the wide range of parameters and the accuracy of spectral measurement with the flexibility of a portable solution for measurements in the field.

- For officially recognised results
- Battery powered
- Lighter than comparable devices
- Robust and waterproof




## Reliable and flexible in all applications

- Wavelength range: 340–800 nm
- For all common cuvette sizes
- For wastewater, drinking water, surface water or industrial water and quality controls

## Robust and practical for all situations

The portable DR1900 combines the highest level of accuracy for everyday use. The housing, made from shock-proof and impact-resistant plastic, offers a non slip surface. The buttons are large and robust, so every touch hits its target.



*With this practical backpack, you can comfortably carry the DR 1900 with you wherever you go and keep your hands free.*

DOC032.52.20119.Aug15

**HACH**<sup>®</sup>  
Be Right™